



TEMPORARY MATERIALS

01_ FILLING

02_ CEMENT

03_ CROWN

FROM TEMPORARY FILLING TO UPPER STRUCTURE HOLE FILLER

01

FOR DENTAL IMPLANT ALL IN ONE SYSTEM



SIMPLE SPEEDY SAFE

Easy to handle
and to remove
in one piece

Time-saving!
Do not using
instrument for
making occlusion

Excellent
sealing
ability

Once-Fil Flow with a great flowability is easily applied to patients.

The elastic force and strength is great for sealing the margin.

Available to remove at one time and provide the strong masticatory forces.

SEALING TEST



Once-Fil Flow

FLOW



CONTENTS

- #101-1 Once-Fil Flow Blue 2g x 5syringes
- #101-2 Once-Fil Flow Yellow 2g x 5syringes
- #101-3 Once-Fil Abutment 2g x 5syringes

Once-Fil HV

HV



CONTENTS

- #102-1 Once-Fil HV Blue 3g x 3syringes
- #102-2 Once-Fil HV Yellow 3g x 3syringes

02

STRONG RETENTION EASY REMOVAL

HIGH RETENTION
EASY CLEAN-UP
EASY HANDLING

PRODUCT	WORKING TIME	SETTING TIME
ANY-TEMP NE	1'00''	3'00''
ANY-TEMP CLEAR	1'30''	5'00''


Advantages

- Long-term cementation
- Easy to remove the excess cement
- Easy to use with Automix type
- Excellent radiopaque
- Low thickness

Application



How to use



Any-Temp Clear



For long term temporary

CONTENTS

- #202-1 Any-Temp Clear 6g x 1syringe
- #202-2 Automixing Cement Tips

Any-Temp NE



For temporary

CONTENTS

- #201-1 Any-Temp NE 10g x 2syringes
- #201-2 Any-Temp NE 10g x 3syringes


STILL USE SMELLY MATERIAL?

Ezi-Crown®

Light cured, temporary crown and bridge resin




Single Crown



Case by Dr. HoYeol Jang
HoYeol Jang Dental Clinic

Bridge Case (Duplication tech)



Case by Dr. HoYeol Jang
HoYeol Jang Dental Clinic

TEMPORARY CROWN COMPLETED IN LESS THAN 5 MIN.



NO SMELL



SAVE TIME



SAVE COST

Application

1. TEMPORARY CROWN & BRIDGE
2. TEMPORARY PROsthESIS RELINING
3. PROsthESIS DUPLICATION
4. IMMEDIATE IMPLANT RESTORATION


Relining Case



Contents :

- #701 EZI-CROWN START KIT
EZI-CROWN A2 15g x 2EA
EZI-CROWN A3 15g x 2EA
EZI-BLOCK 1.2ml x 2syringes / Disposable tips
EZI-SEP 5mL x 1syringe
EZI-CLEAR 50ml x 1EA
- #702 EZI-CROWN REFILL 15g x 3EA
- #703 EZI-CROWN SINGLE 15g x 1EA

Immediate implant restoration



Case by Dr. HoYeol Jang
HoYeol Jang Dental Clinic