

	One-Fil	One-Fil Putty	One-Fil Putty Injectable
Usage	Root Canal Sealer	Pulp capping Apicoectomy Replantation Root-end filling	Pulpotomy Perforation repair Apexification
pH	13	12	12
Setting time	< 3.5 hours	< 5min. (Initial)	< 5min. (Initial)
Flowability [mm]	26	-	-
Solubility [%]	< 1	< 1	< 1
Radiopacity [mmAl]	9	10	10

VIDEO 

One-Fil

One-Fil
One-Fil Putty
One-Fil Putty Injectable

Clinical Application

RESEARCH [SCI]



One-Fil



One-Fil Putty



Hi-Bond Universal



www.mediclus.co.kr

MEDICLUS CO., LTD.

Seoul Office	+82-10-8535-2877
--------------	------------------

No. 206, 464, Cheongpa-ro, Jung-gu, Seoul 04510, Korea
Fax : +82-2-6499-1006 / E-mail : sales@mdclus.com

Head Office	Research & development center
-------------	-------------------------------

No. 1210, 134, Gongdan-ro, Heungdeok-gu, Cheongju-si,
Chungcheongbuk-do 28576, Korea

Calcium silicate based Bioceramic materials



