

Any-Etch[®]

Any-Etch[®] Classic

37% Phosphoric acid gel in two viscosities and two colors

Strong Bonds
Stable Results



Any-Etch®

Any-Etch® Classic

Targeted viscosity

The viscosity and flow characteristics of phosphoric acid gels can significantly impact clinical efficiency. Some products may exhibit reduced viscosity when passing through the syringe tip, leading to imprecise application, or the material may flow uncontrollably into the cavity, causing inconvenience. Any-Etch® and Any-Etch® Classic are specially formulated with consistent thixotropic viscosity and controlled flow, enabling rapid and effective application tailored to intraoral conditions and treatment strategies.

Targeted color

Any-Etch® (Vivid blue) offers high contrast with tooth structure, ensuring precise application and allowing clear differentiation between etched and non-etched areas before and after rinsing. This helps prevent over-etching and ensures optimal bonding surface preparation. Any-Etch® Classic (Soft blue) flows easily into complex cavity morphologies and tooth structures, spreading smoothly and can be removed with light rinsing. Its gentle gel consistency is ideal for intricate preparations, providing effective etching with minimal effort.

Color



Any-Etch®



Any-Etch® Classic



An artwork painted with Any-Etch®

Indications

For pretreatment of teeth and restorations in adhesive procedures: Composite resin restoration, Adhesive bonding of all-ceramic and composite resin indirect restorations, Splint, Orthodontic appliances, Tooth jewelry, Sealants

Delivery Forms



#304-1 Any-Etch®

Any-Etch 3ml x 5ea



#305-1 Any-Etch® Classic

Any-Etch Classic 3ml x 5ea

MEDICLUS

Headquarter | No.1210, 134, Gongdan -ro, Heungdeok-gu, Cheongju-si, Chungcheongbuk-do, 28576, Korea | +82 43 211 2877

Sales Office | No.206 Brown Stone, 464, Cheongpa-ro, Jung-gu, Seoul, 04510, Korea | +82 10 8385 2877

www.mediclus.co.kr | sales@mdclus.com

