


Canal-Seek™

ROOT CANAL ORIFICE & CRACK INDICATOR


Advantages

- Easy application with a syringe type
- Making root canal orifice visible & crack indicator in a short time
- Accurate cavity indication to prevent secondary caries
- Saving treatment time


Application







1. Apply Canal-Seek to the treatment area



2. After 5~10 sec. rinse with water and dry



3. Treat according to the points intensively dyed in blue with the naked eyes



ROOT CANAL ORIFICE & CRACK INDICATOR

CONTENTS

507 Canal-Seek 3ml x 1ea

Endo@Prep™ EDTA Cream

FDA

FILE CREAM



EDTA CREAM TYPE FDA

CONTENTS

502-1 Endo@Prep EDTA Cream 9g x 2ea

502-2 Endo@Prep EDTA Cream 5g x 1ea