

Any-Com®

FDA

NANO-HYBRID COMPOSITE

Light-cured, radiopaque nano-hybrid composite for the restorative treatment of anterior and posterior teeth.


Applications

- Direct fillings of class I, II, V of posterior teeth
- Direct fillings of class III, IV, V of anterior teeth
- Repairs of veneered parts of fixed prosthetic restorations

Advantages

- Excellent radiopacity
- Minimize polymerization shrinkage
- Excellent handling without stickiness and optimum viscosity
- Depth of cure : 2.5mm (10 seconds)
- Variety of shade choices
- Containing zirconia nano fillers
- Improved abrasion resistance by containing zirconia nano filler

Physical properties



SHADE

A1	A2	A3	A3.5	A4	B1	B2	C2	D2	OA2	OA3	TL	WT
									OP	A2E	A3E	



NANO-HYBRID COMPOSITE

CONTENTS

# 301H	Any-Com Universal Kit
	Any-Com 4g x 6ea
	Any-Etch 3ml x 1ea, Hi-Bond Universal 5ml x 1ea
# 301-00	Any-Com Refill 4g x 1ea
# 302H	Any-Com Kit
	Any-Com 4g x 5ea
	Any-Etch 3ml x 1ea, Hi-Bond Universal 5ml x 1ea