

Fast and simple gingival retraction and hemostasis with aluminum chloride

Applications

- Gingival retraction and hemostasis prior to impression taking
- Bleeding control during crown and bridge preparations
- Moisture and bleeding management before bonding or cementation
- Gingival management for Class II and V restorations
- Temporary gingival displacement for improved clinical visibility

Advantages

- Less time, better 2 in 1 effectiveness
- Rapid and reliable bleeding control through the astringent action of aluminum chloride
- Uniform paste application ensures excellent mucosal adhesion and dispersion
- Comfort for both clinicians and patients
- Reduced discomfort, fast and clean rinsing
- Two viscosities and shades to support treatment strategies

Results

Clinical Case by Dr. Min-Sik Choi, Seoul Smart Dental Clinic



Before Treatment



Apply One-Stop™ & Dry



After Treatment

Medium Flow (Yellow)

- For shallow to medium gingival sulcus
- Expansion and drying of the sulcus prior to digital or conventional impression taking
- Suitable for temporary bleeding control

Low Flow (Blue)

- For deep sulcus or thick gingival tissue
- Recommended as a substitute for the second cord in the double-cord technique

Lorem ipsum



GINGIVAL RETRACTION & HEMOSTASIS PASTE

CONTENTS

# 1101-1	One-Stop Blue 0.85g x 7ea
# 1101-2	One-Stop Blue 0.85g x 2ea
# 1102-1	One-Stop Yellow 0.85g x 7ea
# 1102-2	One-Stop Yellow 0.85g x 2ea