

Any-Flex™ (Heavy, Light, Bite, Clear)

IMPRESSION MATERIAL

Any-Flex is hydrophilic vinyl polysiloxane impression materials in the types of Heavy, Light, Bite, Clear.

Heavy body

- It is used for one-step taking of whole teeth or individual tooth with light body
- It is used with light body to make a precise impression

Light body

- It is used for one-step (light body and heavy body) or two-step (light body and putty) impression taking of individual or whole teeth in the mouth

Bite

- It is the bite registration impression material to measure the occlusal surface.
- It is a short polymerization time and a high final hardness
- It is easy to use and remove the low elastic deformation which is suitable for bite registration

Clear

- It is used for one-step taking of several teeth
- Own characteristic allows penetration of dental curing light's source, therefore it is possible to cure the inside material

Time of application

Items	Working Time 23°C	Setting Time (Intra-oral, 23°C)	Shore A
Heavy body	1:00	3:30	61
Light body	1:00	3:30	46
Bite	0:20	0:45	88
Clear	0:25	0:60	60



Any-Flex™ (Heavy, Light, Bite, Clear)

CONTENTS

# 901-1 Any-Flex Heavy Body 50ml x 4ea	# 902-1 Any-Flex Light Body 50ml x 4ea
# 901-2 Any-Flex Heavy Body 50ml x 10ea	# 903-2 Any-Flex Bite Blue 50ml x 2ea
	# 904-1 Any-Flex Clear 50ml x 4ea