

One-Fil® (원-필)

FDA

칼슘실리케이트 계열 바이오세라믹 실러

우수한 밀폐력으로 엔도솔루션을 제공하는 칼슘실리케이트 근관충전용 실러

특징

- 우수한 흐름성으로 편리하게 사용가능
- 우수한 체적안정성과 밀폐력
- 방사선불투과성 900% 이상
- 우수한 생체친화성 및 항균성
- Canal wall과의 탁월한 적심성

경화시간

- 30min in 100% humidity condition. (ISO6876:2012)
- 3.5 hours in normal root canal.

pH



결과



치료 전



치료 후



15개월 후



3개월 후



3년 후



4년 후



Aust Endod J 2021

ORIGINAL RESEARCH

Physicochemical properties and cytocompatibility of newly developed calcium silicate-based sealers

Min-Gyu Park, DDS¹; In-Ryoung Kim, PhD^{2,3}; Hyung Joon Kim, PhD^{2,4}; Sang Won Kwak, DDS, MS, PhD^{1,5}; and Hyeon-Cheol Kim, DDS, MS, PhD^{1,3}



¹ Department of Conservative Dentistry, Dental Research Institute, Pusan National University School of Dentistry, Yangsan, Korea

² Department of Oral Anatomy, Pusan National University School of Dentistry, Yangsan, Korea

³ Dental and Life Science Institute, Pusan National University, Yangsan, Korea

⁴ Department of Oral Physiology, Pusan National University School of Dentistry, Yangsan, Korea

⁵ Periodontal Diseases Signalling Network Research Center, Pusan National University, Yangsan, Korea



Cell Attachment Evaluation (One-Fil)



칼슘실리케이트 근관충전용 실러 FDA

CONTENTS

508-1 One-Fil 2g x 1ea

508-2 One-Fil 0.5g x 1ea

소비자가격

66,000원

22,000원